

**NEW**



**LYMPHEDEMA  
DEPOT LTD**

**LR** **SOARIS**  
COLLECTION  
People.Health.Care

## TributeWrap™

Adjustable, variable-compression foam garments for managing edema and fibrosis.  
Standard, "off-the-shelf" sizing for hands, arms and lower leg.



Wrist to Axilla



Glove



Below Knee

# TributeWrap™

Effective, easy-to-use, and comfortable

Adjustable, variable-compression foam garments to help manage edema, lymphedema and fibrosis.

TributeWrap™ is an off-the-shelf garment that comfortably adjusts to the users' limb and lifestyle. This versatile garment can be used during the evening, night, or low activity — such as working at a desk job or resting on the couch.

While the user is at rest, its proven foam technology works hard to break up hardened or fibrotic tissue to help support soft and healthy skin care. Most importantly, it allows the wearer to maintain the progress made during therapy, while also allowing the flexibility needed to support further reduction.

- **Durable straps** help users to adjust their garment for the right fit
- **Donning tabs** make it easy to get a better grip and pull for improved contour to the limb
- **Powered by Coolcore®** to rapidly wick away sweat and moisture and regulate body temperature
- **Chevron channeling** is designed based on MLD principles of directing fluid to collateral pathways
- **Low profile height** helps provide a comfortable garment experience

## TributeWrap™ Glove

Size	Circumference (cm)		Orientation	Item No.
	Palm (B)	Wrist (C)		
Small	16–20	16.5–21	Left	
			Right	
Medium	18–22	17.5–22	Left	
			Right	
Large	20–24.5	18.5–23	Left	
			Right	

## TributeWrap™ Wrist to Axilla

Size	Circumference (cm)			Arm Length (C–G)	Orientation	Item No.
	Wrist (C)	Elbow (E)	Axilla (G)			
Small	16–21	22–28	25–31	43–47	Left	
					Right	
Medium	17–22	26–32	30–37		Left	
					Right	
Large	18–23	30–36	35–44		Left	
					Right	

## TributeWrap™ Below Knee

Size	Circumference (cm)		Length (cm)		Orientation	Item No.		
	Ankle (B)	Calf (D)	Foot (A–Y*)	Calf (Y*–D)				
Small	24–30	30–36	19–22	40–44	Left			
					Right			
Medium	27–33	33–41.5			19–22	40–44	Left	
							Right	
Large	30–36	36–45			19–22	40–44	Left	
							Right	

\*Y landmark is at the base of the heel.