

# Comparing Solaris Collection TributeNight™ and Caresia™ Garments for the Lower Extremities

## TributeNight™ Full Leg, Half Leg, Foot, Chap, Pants, Shorts



## Caresia™ Below Knee, Thigh and Foot



➤ Custom measured, designed and manufactured for each individual patient.

➤ Colour Options: Black, Slate Grey, Purple, or Raspberry Pink

➤ Bi-directional gradient compression.

➤ Provides approx. 20-30mmHg.  
Addition of optional Outer Jacket adds approx. 10mmHg for 30-40mmHg.  
Variable Compression Jackets are also available.

➤ Contains mixed density, medical grade foam particles ("chip" foam).

➤ Quilted chevron shapes encourage fluid flow laterally and proximally, routing fluid away from dysfunctional areas. Some styles also available in vertical channels option.

➤ Foam particles nestle into skin to provide local tissue stretch, stimulate lymphatics and help soften protein-thickened tissue and resolve fibrosis.

➤ Can begin at toes or ankle; end at knee, mid-thigh, groin or waist. Custom modifications such as zippers, Velcro® closures, pull-up loops and non-skid pads are available. Also available in bilateral, custom-length pants and shorts styles.

➤ Sure-Fit Guarantee™ assures initial fit, covers any material or manufacturing defect, and includes one complimentary reduction alteration (1-Year)

➤ Individually custom manufactured; delivery typically takes 8-12 business days; \$10 flat-rate for insured, expedited Shipping & Handing to any Canadian address

➤ Commonly covered by private insurance plans with a physician's prescription, either under compression garments coverage or durable medical equipment. Zero-rated for GST/HST with a prescription.

➤ Available in a range of standard-sizes for Below Knee, Thigh, and Foot & Ankle.

➤ Comes in Caresia Blue only

➤ Bi-directional gradient compression.

➤ Compression depends on initial fit within range for standard size selected, and can then be customized and adjusted within therapeutic ranges (apx. 20-40mmHg) by external application of short-stretch bandaging.

➤ Contains medical grade foam particles ("chip" foam).

➤ Quilted chevron shapes encourage fluid flow laterally and proximally, routing fluid away from dysfunctional areas.

➤ Foam particles nestle into skin to provide local tissue stretch, stimulate lymphatics and help soften protein-thickened tissue and resolve fibrosis.

➤ Full leg coverage in two pieces: toe-to-knee plus knee-to-groin; Below Knee and Thigh garments are available separately and can be worn either individually or together. Also available in a foot & ankle "boot" style.

➤ Guaranteed to be free of manufacturing defects. All sales final, but sizing exchanges may be accepted; some conditions apply.

➤ In-stock and ships within 24-hours; delivery typically takes 1-5 business days and flat-rate \$10 shipping on any sized order to any Canadian address

➤ Sometimes covered by private insurance plans with a physician's prescription, either under compression garments coverage or durable medical equipment. Subject to applicable Canadian sales taxes.