



# PHYSICIAN'S PRESCRIPTION

Patient: \_\_\_\_\_ Date: \_\_\_\_\_

Physician: \_\_\_\_\_

Signature: \_\_\_\_\_

Prescribing physician contact information or stamp:

## Diagnosis

- Lymphedema (chronic)
- Swelling (acute)
- Chronic Venous Insufficiency
- Combined Lymphedema
- Lipedema

## Area Affected

- Right Arm     Left Arm
- Right Leg     Left Leg
- Right Foot     Left Foot
- Other \_\_\_\_\_

## Treatment

- Solaris Tribute (30mmHg) → →  + Outer Jacket (+10mmHg)
- Solaris ReadyWrap                     + Variable Compression Jacket
- Solaris Caresia                             + Modifications as necessary

Comments:



# PHYSICIAN'S SOLARIS PRODUCT GUIDE

## TRIBUTE

The Solaris Tribute is a custom-made therapeutic garment that is designed to be worn at night as a more consistent and effective alternative to nighttime bandaging. The quilted, directional channels sewn into the fabric are filled with medical-grade foam chips which provide tissue stretch, stimulating the lymphatics and breaking down fibrosis. The channels guide fluid to areas rich in potential for the formation of collateral pathways. The material is breathable, lightweight, and thermo-releasing, with anti-microbial properties to keep the skin cool, dry and clean. Tribute can be designed and created for any part of the body.

**Compression:** Apx. 30mmHg distally decreasing proximally to apx. 20mmHg.

**Indications:** Chronic swelling of any kind. If wounds are present they must be dressed.

**Contraindications:** Include but are not limited to: arterial insufficiency, untreated CHF, obstructive edema, deep vein thrombosis, active infection, absence of sensation, any condition where increased fluid return is a concern.

## TRIBUTE - OUTER JACKETS & VARIABLE COMPRESSION JACKETS

**Outer Jackets** are simple, mildly compressive sleeves to be worn over the Tribute garment.

**Compression:** Adds apx. 10mmHg of compression when worn with Tribute.

**Variable Compression Jackets** are custom-made of strong, low-stretch, Velcro®-receptive material with an overlapping, adjustable tabbed closure system. Compression is set in consultation with the therapist to provide customized compression and accommodate girth fluctuations, particularly reductions.

## READYWRAP

ReadyWrap is a medical compression binder of low elasticity. It is a valuable tool in cases where swelling is aggressive or the venous system requires support. Properly sized and applied it achieves therapeutic ranges of compression for venous support. Its design provides a 50% overlap of straps, which is recommended in clinical bandaging techniques for graduated compression. Currently available for lower extremity applications only (foot, lower leg) where standing pressures exacerbate existing conditions.

**Compression:** Variable compression; therapeutic ranges (20-60 mmHg)

**Indications:** Chronic swelling of any kind. If wounds are present they must be dressed.

**Contraindications:** Same as for Tribute, above.

## CAREZIA COMPRESSION BANDAGING LINER GARMENTS

Caresia's construction features quilted, directional channels filled with medical grade foam chips as found in the Tribute. Caresia is available in a range of standard styles and sizes, and compression is applied and customized by the patient and/or their therapist with an overwrap of short-stretch bandaging. Caresia is extremely useful in both clinical treatment and home-care phases of lymphedema therapy and for addressing any acute edema, particularly where fibrosis is present.

**Compression:** Variable compression; therapeutic ranges (20-50 mmHg)

**Indications:** Chronic swelling of any kind. If wounds are present they must be dressed.

**Contraindications:** Same as for Tribute, above.