



## Solaris Product Guide for Venous Insufficiency

### Solaris Caresia Below-Knee Bandage Liner LE-AD



Soft Caresia Bandage Liners, featuring the same foam chip-filled quilted directional channels as Tribute, provide comfortable and easy-to-apply compression for daytime or overnight wear, making symptom management easier and more effective.

### Solaris Tribute Below-Knee Custom LE-AD-V



For patients who need a more precise fit, the custom made Solaris Tribute provides perfection. Vertical Channels from toes to knee provide directional flow. This style is appropriate for clients who have intact lymph nodes, however, exhibit lymphatic insufficiency, edema, or indurated tissue below the knee. The vertical channels provide excellent ankle definition, effectively addressing ankle swelling.

Learn more at our website:  
[www.LymphedemaDepot.com](http://www.LymphedemaDepot.com)

### Solaris ReadyWrap and Solaris Silver Liners



Chronic venous insufficiency responds well to compression. ReadyWrap is a low-stretch, adjustable compression therapy garment. ReadyWrap employs the 50% overlap pattern of traditional short-stretch bandaging, but its Velcro™ closure system is much quicker and easier for the patient to don and doff. ReadyWrap is available in foot, calf, knee and thigh units that can be worn separately or combined for full leg coverage.

Every ReadyWrap garment is constructed of comfortable, breathable and moisture-wicking Breathoprene™ fabric.

Solaris Silver Liners are manufactured with a blend of materials that result in a soft, comfortable and breathable fabric with anti-microbial, anti-bacterial and anti-static properties due to the presence of silver fibres.

Wearing Solaris Silver Liners not only offers the wearer comfort and protection but also protects compression garments from direct contact with the skin, helping to extend the life of garments and keeping them fresher longer.

